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

9,217 Estimated number of cell
5,509 Median UMI counts per cell
2,057 Median genes per cell
46,096 Mean reads per cell

人PBMC

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
10,528 Estimated number of cell
3,280 Median UMI counts per cell
1,393 Median genes per cell
36,807 Mean reads per cell

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
10,163 Estimated number of cell
5,585 Median UMI counts per cell
2,475 Median genes per cell
33,308 Mean reads per cell

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
10,387 Estimated number of cell
4,226 Median UMI counts per cell
1,498 Median genes per cell
18,218 Mean reads per cell

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
9,045 Estimated number of cell
4,775 Median UMI counts per cell
2,007 Median genes per cell
31,428 Mean reads per cell

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
12,979 Estimated number of cell
2,165 Median UMI counts per cell
1,141 Median genes per cell
17,655 Mean reads per cell

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
8,119 Estimated number of cell
2,692 Median UMI counts per cell
1,026 Median genes per cell
32,567 Mean reads per cell

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7,822 Estimated number of cell
4,773 Median UMI counts per cell
1,785 Median genes per cell
28,984 Mean reads per cell

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
- [1] Mazid MA. et al. Rolling back human pluripotent stem cells to an eight-cell embryo-like stage. *Nature*. (2022)
- [2] Han L. et al. Cell transcriptomic atlas of the non-human primate *Macaca fascicularis*. *Nature*. (2022)
- [3] Zhu, L. et al. Single-cell sequencing of peripheral mononuclear cells reveals distinct immune response landscapes of COVID-19 and influenza patients. *Immunity*. (2020)
- [4] Zhao, Z. et al. Evolutionarily conservative and non-conservative regulatory networks during primate interneuron development revealed by single-cell RNA and ATAC sequencing. *Cell research*. (2022)
- [5] Li, Q. et al. A single-cell transcriptomic atlas tracking the neural basis of division of labour in an ant superorganism. *Nature Ecology & Evolution*. (2022)
- [6] Liu, K. et al. Bilineage embryo-like structure from EPS cells can produce live mice with tetraploid trophectoderm. *Protein & Cell*. (2022)
- [7] Sun, S. et al. *Znht1* controls meiotic initiation in male germ cells by coordinating with *Stra8* to activate meiotic gene expression. *Developmental Cell*. (2022)
- [8] Liu, C. et al. Spatiotemporal mapping of gene expression landscapes and developmental trajectories during zebrafish embryogenesis. *Developmental Cell*. (2022)
- [9] Ding, R. et al. Single-cell transcriptome analysis of the heterogeneous effects of differential expression of tumor PD-L1 on responding TCR-T cells. *Theranostics*. (2022)
- [10] Wang, W. et al. Multiregion single - cell sequencing reveals the transcriptional landscape of the immune microenvironment of colorectal cancer. *Clinical and translational medicine*. (2021)
- [11] Li, Z. et al. Molecular mechanisms governing circulating immune cell heterogeneity across different species revealed by single - cell sequencing. *Clinical and Translational Medicine*. (2022)
- [12] Chen, D. et al. Single - cell atlas of peripheral blood mononuclear cells from pregnant women. *Clinical and translational medicine*. (2022)
- [13] Ding, R. et al. Single-cell transcriptome analysis of the heterogeneous effects of differential expression of tumor PD-L1 on responding TCR-T cells. *Theranostics*. (2022)
- [14] Li, Y. et al. *Npas3* deficiency impairs cortical astrogenesis and induces autistic-like behaviors. *Cell Reports*. (2022)
- [15] Xu, X. et al. Single Nuclear RNA Sequencing Highlights Intra-Tumoral Heterogeneity and Tumor Microenvironment Complexity in Testicular Embryonic Rhabdomyosarcoma. *Journal of Inflammation Research*. (2022)
- [16] Liu, X. et al. Single-Cell Atlas of the Chinese Tongue Sole (*Cynoglossus semilaevis*) Ovary Reveals Transcriptional Programs of Oogenesis in Fish. (2022)
- [17] Shang, Z. et al. Integrative single-cell RNA-seq and ATAC-seq analysis of mesenchymal stem/stromal cells derived from human placenta. *Frontiers in cell and developmental biology*. (2022)

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